

LIGHT WHERE YOU NEED IT



14
KLX AT 4"



3 HR.
RUN TIME



IPX5
RESISTANT



**MADE IN
THE USA**

STERILE EO
STERILIZED USING
ETHYLENE OXIDE



**NOT MADE
WITH
NATURAL
RUBBER
LATEX**



5 PACK
BOXES



**SINGLE USE
DISPOSABLE**

Scintillant[®] Suction



Part No. 2658-01-0002

The Scintillant[®] Suction is a handheld surgical light assembled to a Yankauer style suction tip designed to combine fluid aspiration with bright LED illumination. The Suction is ready for use directly from the package, enabling users to immediately navigate obstructions safely, clearly visualize anatomy, and aspirate fluids easily.

CE 0344 3-Year Shelf Life

Scintillant[®] Retractor



Part No. 2658-01-0004

The Scintillant[®] Retractor is a handheld surgical light affixed to a low profile flat retractor designed to combine moderate blunt retraction of tissue with bright LED illumination. The ergonomic design enables users to safely navigate critical obstructions, clearly visualize anatomy and efficiently retract tissue. The Retractor is ready for use immediately upon removal from its protective packaging.

CE 0344 3-Year Shelf Life

Scintillant[®] Zip-Strap



Part No. 2658-01-0001

The Scintillant[®] Zip-Strap is a handheld surgical light featuring rubber straps for simple attachment to surgical instruments. The non-tethered construction enables users to enhance operating room instruments, clearly illuminate the surgical field and safely navigate critical obstructions.

CE 0344 3-Year Shelf Life

Scintillant® SURGICAL LIGHT

SURGICAL PROCEDURE	Suction	Retractor	Zip-Strap
GYNECOLOGICAL	2658-01-0002	2658-01-0004	2658-01-0001
Abdominal Hysterectomy	●		●
Abdominal Hysterectomy-Obese Pt	●	●	●
Pelvic Floor Reconstruction	●	●	●
Prolapse Procedures	●	●	●
Vaginal Hysterectomy	●	●	●
Vaginal Vault Suspension	●	●	●
UROLOGIC			
Bladder Suspension	●	●	
Cystectomy	●		●
Nephrectomy	●		●
Prostatectomy	●	●	●
Renal Decapsulation			●
COLORECTAL			
Abdominal Perineal Resection	●		●
Lower Anterior Resection	●	●	●
Sigmoid Colon Resections	●	●	●
GENERAL			
Hepatic Procedures	●	●	●
Lymph Node Dissection	●	●	●
Oncological Procedures	●	●	●
ORTHOPEDIC			
Hip Revision	●		●
Total Hip Replacement	●		●
PLASTIC/RECONSTRUCTIVE			
Abdominoplasty	●	●	●
Breast Augmentation	●	●	●
Craniofacial Reconstruction	●	●	●
Reconstruction Mammoplasty	●	●	●
Thoracic	●	●	●
Esophagectomy	●	●	●
Lobectomy	●	●	●
Nipple Sparing Mastectomy	●	●	●
Thoracotomy	●	●	●
SPINE			
Deep Cavity Illumination			●

SEE IFU FOR INTENDED USE AND CONTRAINDICATIONS

2658-80-0044 Rev B